Health & Safety Warnings

* These health & safety warnings are periodically updated for accuracy and completeness. Check www.lenica.ca/warnings for the latest version.



HEALTH & SAFETY WARNINGS: TO REDUCE THE RISK OF PERSONAL INJURY, DISCOMFORT OR PROPERTY DAMAGE. PLEASE ENSURE THAT ALL USERS OF PEAK COGNITION READ THE WARNINGS RELOW CAREFULLY REFORE USING PEAK COGNITION

A WARNING Before Using Peak Cognition:

- Read and follow all setup and operating instructions provided with Peak Cognition.
- Review the hardware and software recommendations for use of Peak Cognition. Risk of discomfort may increase if recommended hardware and software is not used.
- Your Peak Cognition is not designed for use with any ineligible device, accessory and/or software. Use of an ineligible device, accessory and/or software may result in injury to you or others, and may cause performance issues or damage to your system and related services
- · A comfortable virtual reality experience requires an unimpaired sense of motion and balance. Do not use Peak Cognition when you are tired, need sleep, are under the influence of alcohol or drugs, are hung-over, have digestive problems, are under emotional stress oranxiety, orwhen suffering from cold, flu, headaches, migraines, or earaches, as this can increase your susceptibility to adverse symptoms.
- We recommend seeing a doctor before using Peak Cognition if you are pregnant, elderly. have pre-existing binocular vision abnormalities or psychiatric disorders, or suffer from a heart condition or other serious medical condition

A WARNING Seizures:

Some people (about 1 in 4000) may have severe dizziness, seizures, epileptic seizures or blackouts triggered by light flashes or patterns, and this may occur while they are watching TV, playing video games or experiencing virtual reality, even if they have never had a seizure or blackout before or have no history of seizures or epilepsy. Such seizures are more common in children and young people under the age of 20. Anyone who has had a seizure, loss of awareness, or other symptom linked to an epileptic condition should see a doctor before using Peak Cognition.

A WARNING Children:

Peak Cognition should not be used by children under the age of 13, as young children are in a critical period in visual development. Adults should monitor children (age 13 and older) who are using or have used Peak Cognition for any of the symptoms described below, and should limit the time children spend using Peak Cognition and ensure they take breaks during use. Prolonged use should be avoided, as this could negatively impact hand-eye coordination, balance, and multi-tasking ability. Adults should monitor children closely during and after use of Peak Cognition for any decrease in these abilities.

A WARNING General Precautions:

To reduce the risk of injury or discomfort you should always follow these instructions and observe these precautions while using Peak Cognition:

- . Use Only In A Safe Environment: Peak Cognition produces an immersive virtual reality experience that distracts you from and completely blocks your view of your actual surroundings.
 - Always be aware of your surroundings before beginning use and while using Peak Cognition. Use caution to avoid injury.
 - Use of Peak Cognition may cause loss of balance.
 - . Remember that the objects you see in the virtual environment do not exist in the real environment, so don't sit or stand on them or use them for support.
 - For the safest experience, remain seated or stationary when using Peak Cognition.
 - Serious injuries can occur from tripping, running into or striking walls, furniture or
 - other objects, so clear an area for safe use before using Peak Cognition. Take special care to ensure that you are not near other people, objects, stairs. balconies, open doorways, windows, furniture, open flames, ceiling fans or light fixtures or other items that you can bump into or knock down when using—or immediately after using—Peak Cognition.
 - Remove any tripping hazards from the area before using Peak Cognition. Remember that while using Peak Cognition you may be unaware that people and
 - pets may enter your immediate area. Do not handle sharp or otherwise dangerous objects while using Peak
 - Cognition.
 - Never wear Peak Cognition in situations that require attention, such as walking, bicycling, or driving,

- Make sure Peak Cognition is level and secured comfortably on your head, and that you see a single, clear image.
- Ease into the use of Peak Cognition to allow your body to adjust; use for only a few minutes at a time at first, and only increase the amount of time using Peak, Cognition gradually as you grow accustomed to virtual reality. Looking around and using the input device when first entering virtual reality can help you adjust to any small differences between your real-world movements and the resulting virtual reality experience.
- Do not use Peak Cognition while in a moving vehicle such as a car, bus, or train, as this can increase your susceptibility to adverse symptoms.
- Take at least a 10 to 15 minute break every 30 minutes, even if you don't think you need it. Each person is different, so take more frequent and longer breaks if you feel discomfort. You should decide what works best for you.
- Peak Cognition may be equipped with a "pass-through" feature which permits you to temporarily see your surroundings for brief real world interaction. You should always remove Peak Cognition for any situation that requires attention or coordination.
- Listening to sound at high volumes can cause irreparable damage to your hearing. Background noise, as well as continued exposure to high volume levels, can make sounds seem quieter than they actually are. Due to the immersive nature of the virtual reality experience, do not use Peak Cognition with the sound at a high volume so that you can maintain awareness of your surroundings and reduce the risk of hearing damage.

A WARNING Discomfort

- Immediately discontinue use if anyone using Peak Cognition experiences any of the following symptoms; seizures; loss of awareness; eve strain; eve or muscle twitching; involuntary movements; altered, blurred, or double vision or other visual abnormalities: dizziness: disorientation: impaired balance: impaired hand-eve coordination; excessive sweating; increased salivation; nausea: lightheadedness: discomfort or pain in the head or eyes: drowsiness: fatigue; or any symptoms similar to motion sickness.
- Just as with the symptoms people can experience after they disembark a cruise ship. symptoms of virtual reality exposure can persist and become more apparent hours after use. These post-use symptoms can include the symptoms above, as well as excessive drowsiness and decreased ability to multi-task. These symptoms may put you at an increased risk of injury when engaging in normal activities in the real world.

 Do not drive, operate machinery, or engage in other visually or physically demanding activities that have potentially serious consequences (i.e., activities in which









- or other activities that require unimpaired balance and hand-eve coordination (such as playing sports or riding a bicycle, etc.) until you have fully recovered from any symptoms. Do not use Peak Cognition until all symptoms have completely subsided for several
- hours. Make sure you have properly configured Peak Cognition before resuming HISA Be mindful of the type of content that you were using prior to the onset of any symptoms
- because you may be more prone to symptoms based upon the content being used. See a doctor if you have serious and/or persistent symptoms.

A WARNING Repetitive Stress Injury:

Using Peak Cognition can make your muscles, joints or skin hurt. If any part of your body becomes tired or sore while using Peak Cognition, or if you feel symptoms such as tingling, numbness, burning or stiffness, stop and rest for several hours before using Peak Cognition again, If you continue to have any of the above symptoms or other discomfort during or after play, stop using Peak Cognition and see a doctor.

A WARNING Interference with Medical Devices:

Peak Cognition can use hardware that contains magnets and can emit radio waves. These can affect the operation of nearby electronics, including cardiac pacemakers, hearing aids and defibrillators. If you have a pacemaker or other implanted medical device, do not use Peak Cognition without first consulting your doctor or the manufacturer of your medical device. Maintain a safe distance between Peak Cognition and your medical devices, and stop using Peak Cognition if you observe a persistent interference with your medical device.

A WARNING Controller:

- Your Peak Cognition may have come with a third-party controller.
- Consult the manufacturer for health and safety warnings for the controller.
- Always use included wrist straps with controllers to secure the controller to your wrist when in use

A WARNING Batteries:

- CHOKING HAZARD. The controllers are not toys. They contain batteries, which are small parts. Keep away from children under 3.
- DO NOT INGEST ÁNY BATTERY. CHEMICAL BURN HAZARD.
- If a battery is swallowed, it can cause severe internal burns and potential perforation of
- esophagus in just 2 hours and can lead to death. Keep batteries out of reach of children. If you think batteries might have been swallowed or placed inside any part of the body, seek
- medical attention, and have your doctor call the battery ingestion hotline at 1-800-498-8666. . If the battery compartment does not close securely, stop using the product and keep it
- away from children. Risk of fire. Batteries could cause thermal event or leak if installed backwards, disassembled. overcharged, punctured, crushed, mixed with used or other battery types, or directly exposed to fire or high temperature. Do not recharge disposable batteries. Replace all batteries in
- a single device at the same time. Keep in original package until ready to use. Properly dispose of used batteries promotly. . Do not allow conductive materials to touch the battery terminals on the devices. Keep
- batteries dry and avoid water intrusion.
- Do not disassemble, pierce, or attempt to modify batteries.

A WARNING Electrical Shock:

To reduce risk of electric shock

- Do not modify or open any of the components provided.
- Do not use the product if any cable is damaged or any wires are exposed.

If a power adapter is provided:

- Do not expose the power adapter to liquid or moisture.
- . Unplug the power adapter before cleaning, and clean only with a dry cloth.
- Keep the power adapter away from open flames and other heat sources.
- Use only the power adapter provided with approved Peak Cognition hardware.

CAUTION Damaged or Broken Device:

- Do not use your device if any part is broken or damaged.
- . Do not attempt to repair any part of your device yourself. Repairs should only be made by an authorized servicer

A CAUTION Contagious Conditions:

To avoid transferring contagious conditions (like pink eye), do not share Peak Cognition with persons with contagious conditions, infections or diseases, particularly of the eyes. skin or scalp Peak Cognition should be cleaned between each use with skin-friendly antibacterial wipes (particularly the lenses) and dried with a microfiber cloth.

CAUTION Skin Irritation:

Peak Cognition is worn next to your skin and scalp. Stop using Peak Cognition if you notice swelling, itchiness, or other irritation. If symptoms persist, contact a doctor,

NOTICE To avoid damage to Peak Cognition:

- DonotleavePeak Cognitionindirectsunlight. Exposure to directsunlight can damagePeak Cognition.
- Do not shine a laser or external light source through the lenses as it may damage the screen.
- Store components in their storage case when not in use to minimize unintentional damage or environmental exposure.